

# FUME HOOD

## TI-12 SERIES

Tempo make fume hood is a type of local ventilation device that is designed to limit exposure to hazardous or toxic fumes, vapors or dusts. Air drawn from lab space into the fume hood then dispersed via pipe to outside atmosphere. Basically it protect the user from inhaling toxic gases.

Available in following type

- Table top / bench top Fume Hood
- Floor mounted fume Hood

### CONSTRUCTION EXHAUST SYSTEM

- Blower with Anti Coregent Coating.
- Low Noise Blower

### MATERIAL OF CONSTRUCTION

- Double Wall Construction
- Inner Body PP & Outer Body M.S. Powder Coated
- Inner Body Fiber & Outer Body M.S. Powder Coated
- Smooth Black Granite Table

### FEATURES

- Interior baffle4s designed to minimize turbulence & optimize containment
- Easy & Smooth operatable toughened glass sash for full view.
- Provides with 9 feet length fabric flexible pipe.
- Two switch sockets provide
- Provides with two gas value
- Provides with LED lighting
- Thick black granite for working



TABLE TOP MODEL



WITH STORAGE CABINET

## OPTIONAL FEATURES

- Storage Cabinet
- Water value with cup sink arrangement
- PP/FRP/HDFE coated pipe
- Flame proof lighting

## TECHNICAL DATA

Models	Internal Size H x W x D in cm	External size H x W x D in cm	Motor Position
TABLE TOP MODEL			
12FM-A	110 x 60.5 x 75	160 x 92 x 112	Inbuilt Blower Motor
12FM-B	122 x 91.5 x 75	225 x 122 x 112	Inbuilt Blower Motor
WITH STORAGE CABINET			
12FM-C	122 x 122 x 75	295 x 155 x 112	Wall Mounted Blower
12FM-D	122 x 152.5 x 75	295 x 185 x 112	Wall Mounted Blower

## APPLICATION

- Chemical Lab
- Food Manufacturer
- Drugs Manufacturer Company
- Paint Manufacturing Company
- Rubber Industry
- Polymer Industry
- Cement Industry
- R & D Labs
- Jewellery Manufacturers

**MUMBAI :** @ temposales@shandilyagroup.com  
**(H.O)** @ sales@tempo.net.in  
☎ +91 9820464003;  
+91 9619741622 ;  
+91 9820840572;  
+91 9152097537



## HYDERABAD

@ Email Id:- tempohyd@shandilyagroup.com  
☎ +91 9959531765

### TEMPO INSTRUMENTS PVT. LTD.

TEMPO HOUSE  
Top syringe compound, W.E. Highway, Pandurang wadi, behind  
Samrat Hotel, Mira Road (East), Thane-401104  
Tel : +91- 9820464003 ; Email : sales@tempo.net.in

### Head Office

10-11, Prospect Chamber Annexe, 317-21, Dr. D.N. Road, Fort,  
Mumbai-400001 INDIA  
Tel : +91-022-65206102-03, 65208849-50 ; Fax: +91-022-22045352