

LEAK-TEST APPARATUS

TI - 141 SERIES

Fully Automatic:

Tempo make fully automatic leak test apparatus. It is used to check the quality of packing and quality control process.

Fully automatic leak test apparatus gives Reliable, Economical and fully calibrated solution in quality assurance for testing of tablet strip, sachet, candies etc.

Features:

- LCD digital display.
- Soft touch keypad.
- Inbuilt printer interface for Dot Matrix Printer.
- Programmable vacuum hold time upto 99 mins.
- Automatic vacuum release.
- Made up of SS304
- Polycarbonate vacuum desiccator.
- Very low noise model.



Technical Data :

Model No.	141 A	141B
Dessicator Diameter	150mm	300mm
Rating	230V	230V
Vacuum Level	550mmhg	550mmhg
Timer	Inbuilt Digital Timer	Inbuilt Digital Timer
Dessicator	Polycarbonate	Polycarbonate
Vacuum Display	Digital Vacuum Display	Digital Vacuum Display
Vacuum pump	Oil free Diaphragm pump	Oil free Diaphragm pump
Body	S.S. 304	S.S. 304
Print	Printer port	Printer port

Application

- Drugs Manufacturer (Blisters and small sachets containing tablets)
- Paint Manufacturing Company
- Food Manufacturer (Sugar, Spices, Granulates).
- Bislery water Bottle Manufacturing

MUMBAI : @ temposales@shandilyagroup.com
(H.O)

@ sales@tempo.net.in

+91 9820464003;
+91 9619741622 ;
+91 9820840572;
+91 9152097537



ISO 9001:2015



ISO 9001
REGISTERED
ISO 13485:2012

HYDERABAD

@ Email Id:- tempohyd@shandilyagroup.com

+91 9959531765

TEMPO INSTRUMENTS PVT. LTD.

TEMPO HOUSE

Top syringe compound, W.E. Highway, Pandurang wadi, behind

Samrat Hotel, Mira Road (East), Thane-401104

Tel : +91-9820464003 ; Email : enquiry@tempo.net.in / sales@tempo.net.in

Tempo

www.tempoinstruments.com