

Tempo

Controlled Temperature Solution

MUFFLE FURNACE

1200°/1000°C

TECHNICAL DETAILS

TI - 58 SERIES : (Maximum Temperature 930° C)

Technical Specification :

Models	Internal Size D x W x H (cm)	Rating (KW)	Max Continuous temp.
TI-58A	23 x 10 x 10	1.6	900°C
TI-58B	30 x 15 x 15	3.3	900°C



TEMPO INSTRUMENTS PVT. LTD.

MUMBAI (H.O):

Top Syringe Compound, W. E. Highway, Pandurang Wadi,
Mira Road (East), Thane 401 104. Maharashtra, INDIA.

@ sales@tempo.net.in; enquiry@tempo.net.in
temposales@shandilyagroup.com



ISO 9001:2015



ISO 13485:2012



HYDERABAD

@ Email Id:- tempohyd@shandilyagroup.com

@ +91 9959531765 / 7893995434