

1200°C MUFFLE FURNACE WITH SHELL AND THYRISTOR (Continuous Max Temp 1150°C)

Model TI-59 HTA / HTB / HTC / HTD

“TEMPO” make Muffle Furnace are designed & manufactured to comply with requirements of Laboratory & Industrial applications. It is used for engineering tests of soils and aggregates, cement testing, ashing organic and inorganic samples, gravimetric analysis, ignition tests, etc. “TEMPO” make Muffle Furnace are constructed in Horizontal & Rectangular model. Insulation is with high grade insulation material.

Technical Specification :

Models	Internal Size D x W x H (cm)	Rating (KW)	Max Continuous temp.
TI-59HTA	23 x 10 x 10	1.8	1150°C
TI-59HTB	30 x 15 x 15	3.5	1150°C
TI-59HTC	30 x 20 x 20	3.5	1150°C
TI-59HTD	45 x 22.5 x 22.5	7	1000°C

Technical Data :

Max. Temp. : 1200°C (Continuous max 1150°C
& 1000°C for 59HTD)
Supply : 230V, AC



TEMPO INSTRUMENTS PVT. LTD.

MUMBAI (H.O):

Top Syringe Compound, W. E. Highway, Pandurang Wadi,
Mira Road (East), Thane 401 104. Maharashtra, INDIA.

@ sales@tempo.net.in; enquiry@tempo.net.in
temposales@shandilyagroup.com

+91 9820464003;
+91 9820840572;
+91 9152097538



ISO 9001:2015



ISO 13485:2012



HYDERABAD

@ Email Id:- tempohyd@shandilyagroup.com
+91 9959531765 / 7893995434