

STATIC PASS BOX

TI - 80 SERIES

Static Pass Box : Static Pass box Box equipment are made of high grades stainless steel (304 / 316) with international quality standard components including electrically interlocked. This pass box is designed to minimize traffic and contamination entry into the clean room. Enables parts, tools and other work items to be passed to and from the room. This pass box unit acts as an air lock device preventing ambient air from entering, or clean air from exiting the clean room.



Features of Static Pass Box

- High grade stainless steel construction.
- Precision designing and finish.
- Electrical interlocked.
- Floor and wall mounted design.
- Buzzer indication to know the material is kept inside.
- Maximum chemical resistance construction.
- UV light along with hour meter.
- LED Light
- Indicating lamp to know the door lock-unlock condition.
- Toughen glass window.
- LCD display

TECHNICAL DATA

Model No	Internal Size	External Size
	H X W X D in cm	H X W X D in cm
TI – 80 A-SPB	45 X 45 X 45	65 X 57 X 53
TI – 80 B-SPB	60 X 60 X 60	81 X 73 X 69
TI – 80 C-SPB	75 X 75 X 75	95 X 87 X 83
TI – 80 D-SPB	90 X 90 X 90	112 X 104 X 100

MUMBAI : @ temposales@shandilyagroup.com
(H.O) @ sales@tempo.net.in
☎ +91 9820464003;
+91 9619741622 ;
+91 9820840572

HYDERABAD

@ Email Id:- tempohyd@shandilyagroup.com
☎ +91 9959531765

TEMPO INSTRUMENTS PVT. LTD.

TEMPO HOUSE
Top syringe compound, W.E. Highway, Pandurang wadi, behind
Samrat Hotel, Mira Road (East), Thane-401104
Tel: +91- 9820464003 ; Email : sales@tempo.net.in



ISO 9001:2015



ISO 13485:2012

Head Office

10-11, Prospect Chamber Annexe, 317-21, Dr. D.N. Road, Fort,
Mumbai-400001 INDIA
Tel: +91-022-65206102-03, 65208849-50 ; Fax: +91-022-22045352